

SAS • HOLDER & TRIPOD



5215

NORMAL FLOOR TRIPOD

- Extendable from 100 cm up to 280 cm high
- It can be used for air conditioned outlet port



42704



89722



STAINLESS STEEL CLEAN ROOM FLOOR TRIPOD

- Robust tripod
- S/S to avoid particle emission in Clean Room
- To collect bioaerosol at the level of breathing
- Fixed high at 168 cm

COMPACT "FLOOR" TRIPOD

STAINLESS STEEL WALL/SURFACE/CART TRIPOD

- For a better free space inside the Clean Room
- For an easy transfer on board of wheel cart
- Fixed high 22 cm

20914



5214

TABLE TRIPOD

- Compact unit
- To be used on the table

WHAT TO ORDER

CODE	DESCRIPTION	PACKAGING
20914	"SAS Holder" s/s wall/surface/cart tripod	1/box
89722	"Clean-Tripod" s/s Clean Room floor tripod	1/box
5215	"Telesco Floor" normal floor tripod (extendable up to 280 cm high)	1/box
42704	"Compact Digi" compact floor tripod (extendable from 50 cm to 165 cm high) with quick disconnect system	1/box
5214	"Telesco Table" table tripod	1/box