

Sas Super

“SAS SUPER ISOLATOR”



**Critical air monitoring in isolator
and barrier environments**

“SUPERSAS ISOLATOR”

For critical air monitoring in isolator and barrier environments

It consists of a programmable unit which remains outside the controlled area and an independent sampling head which uses **RODAC** contact Plates or Standard 90 mm Petri dishes.



- To be placed anywhere within the isolator or barrier environment
- No disturbance of the unidirectional air flow
- Fast sampling (up to 1000 litres of air in ten minutes or six minutes)
- Fully programmable (air volume, delays, multiple cycles) from outside the isolator
- Compatible with most commonly used sanitizing agents
- Printed documentation (sampling location, operator name, sampling parameters)
- Any type of RODAC contact plates or standard 90 mm Petri dishes may be used.

WHAT TO ORDER

Catalogue Number	Description	Quantity / Pck
43216	«SuperSAS Isolator 100» for 100 lts / minute air flow for RODAC Contact Plates (without aspirating head)	1
15503	“Drilled SS-ST” Stainless / steel aspirating head for 55 mm RODAC Contact Plates	1
19119	“Drilled SS-Petri” Stainless / steel aspirating head for standard 90 mm Petri dishes, complete of adaptor	1
18201	“Revit” Battery Charger 220 Volts	1
19248	“Revit” Battery Charger 110 Volts	1
32879	“Politecnico” validation system	1