

Press Release

FIVE THOUSAND LITRES OF AIR IN FOUR MINUTES

Ref. 81/08 En

Milan, 12/12/2008

“SAS-Cyclone” is the name of the new “International PBI” microbiological air sampler indicated to collect large volumes of air in short times.

Up to 5 cubic metres of air can be aspirated in less than 4 minutes.

The micro-organisms are collected in a sterile fluid for the subsequent analysis by Real Time PCR or chips, if the results should be obtained in a few minutes, or by traditional microbiological testing.

The “SAS-Cyclone” is battery or main operated, it uses sterile disposable collecting container, and it is transferred on the field by a Army certified portable case.

The specific version “SAS-Cyclone ARMY” is dedicated to be used in the battle field for the biological war defence.

