

Press Release

COMPLIANCE WITH GOOD AIR MONITORING PRACTICE

Ref. 85/08 En

Milan, 19/12/2008

The "GAP" (Good Air-Monitoring Practice) requests that all the material used for the microbiological monitoring of the environment is sterile and the inspector authorities ask to document the sterilization compliance.

The "Dispo-Head" sterile disposable aspirating head for microbiological air samplers of the "SAS" sampler family has been produced for this purpose.

Each "Dispo-Head" from International PBI is individually packaged and include a sterilization certificate.

If used with the "DUO-SAS-360", 1 cubic metre of air can be samples in about 3 minutes.



"Dispo-Head" sterile disposable aspirating
"Duo SAS 360" Dual Head portable microbiological air sampler