

Press Release

INDOOR AIR QUALITY MONITORING

Ref. 31/08 En

Milan, 06/05/2008

The microbiological quality of the indoor air is an important parameter to protect the health of patients and staff in hospital, students in schools and universities and in all closed crowded public environments and buildings.

The "**SAS-ISO**" from "International PBI" is the new generation of **microbiological air samplers** with several useful features to obtain correct evaluation of the environmental conditions like interval sampling, eight hours autonomy, up to 300 sampling cycles memorized, air speed sensor, automatic calibration reminder, possibility to use 55 mm, 84 mm RODAC plates, standard 90 mm disposable or re-usable Petri dishes, data transfer to P.C. or printer.

