

Press Release

AIR MONITORING IN CLEAN ROOM AND ISOLATOR

Ref. 41/08 En

Milan, 30/06/2008

All analytical materials that are introduced in Clean Room for microbiological testing should be sterile and, according to the GLP and Control authorities, a document of certification should prove that they are correctly sterilized.

This rule is also valid for the aspiration head that is used with microbiological air samplers. The cycle of air sampling should use a documented sterile aspirating head.

“Dispo-Head” is a disposable aspirating head for microbiological air monitoring with “SAS” sampler, which has been developed and produced by International PBI to simplify the workload of the microbiologist.



“Dispo-Head”, made in recyclable antistatic ABS resin plastic, is triple packed, irradiated and comes with certification of sterilisation, conformity and endotoxin level (LAL test). It has a shelf-life of five years and it is suitable for any type of SAS Super air sampler that uses diameter 55 mm Contact Plates (RODAC).