

Press Release

BAR-CODE IN AIR SAMPLING

Ref. 46/08 En

Milan, 29/07/2008

The application of Bar Code to the microbiological air monitoring procedures can help to save time, to better control the activity of the operator and to achieve a complete traceability of the test.

The "SAS ISO" microbiological air sampler from International PBI can be used in combination with a Bar Code system to facilitate and register the important steps of air monitoring such as: plate with nutrient media, air sampler, operator, location identifications.

All sampling data are collected and printed at the end of the day / shift.



"SAS Super ISO"

