

## Press Release

# ISOLATOR MICROBIAL MONITORING

*Ref. 10/08 En*

**Milan, 14/02/2008**

Monitoring of atmosphere inside isolators is very critical in maintaining high quality standards. There are risk in using portable environmental samplers because they must be introduced each time you sample.

“[SAS-Isolator](#)” from International pbi is a dedicated instrument for this application.

The sampling aspirating chamber of “SAS-Isolator” is compact, manufactured in stainless steel and separated from the control unit. The compact design save space inside the isolator. The connection is only for power and doesn't compromise the isolator integrity.

It is also very easy to install. The simplicity of the system avoids the extraction and control of air, eliminating risk of contamination.

