

Coldstar

LABORATORY VERTICAL FREEZER

Range of temperature:
from **-5 to -40 °C**

Automatic defrosting for the best use of operator's time

High density wall insulation to reduce the energetic consumption

Epoxy coated steel structure with rounded internal corner for easy cleaning and disinfection

Insulated hinged lid with removable magnetic rubber gasket and heated frame to guarantee a hermetic closure



Control Panel, microprocessor operating for all functions, acoustic and visual LED display

Forced air to obtain the best temperature uniformity and stability inside the cabinet

R404a and HCFC-free refrigeration system

Inside light activated at door opening

Compact sizes for better use of the laboratory space

Two different control panels RED or SU-X

WHAT TO ORDER

Pannel	Code		Temperature Freezer	Capacity (litres)	Dimensions	Weight (kg)	N° of doors	N° of shelves
	RED	SU-X						
	89735	89736	-20 -40 °C	270	600x650x1640 h mm	80	1	2
	89731	89732	-5 -20 °C	370	600x600x1920 h mm	135	1	3
	89737	89738	-20 -40 °C	530	720x700x2000 h mm	150	1	3
	89378	89733	-5 -25 °C	625	720x800x2000 h mm	160	1	3
	89502	89734	-5 -25 °C	1365	1440x800x2000 h mm	230	2	6