

Coldstar Chest

LABORATORY HORIZONTAL FREEZER

Range of temperature: **from -20 to -45 °C**

Manual defrosting

High density wall insulation to reduce the energy consumption

Epoxy coated steel structure with rounded internal corner for easy cleaning and disinfection



Insulated hinged lid with removable magnetic rubber gasket to guarantee a hermetic closure

Control Panel, microprocessor operating for all functions, acoustic and visual LED display

R404a and HCFC-free refrigeration system

Compact sizes for better use of the laboratory space

Two different control panels RED or SU-X

WHAT TO ORDER

Pannel	Code		Capacity (litres)	Dimensions	Weight (kg)
	RED	SU-X			
	89739	89740	60	400x680x930 h mm	60
	89741	89742	320	1280x680x930 h mm	90