

# Superfreeze

## LABORATORY VERTICAL FREEZER

Range of temperature:  
from -10 to -25 °C

High density wall insulation to  
reduce the energy consumption

Epoxy coated steel structure with  
rounded internal corner for easy  
cleaning and disinfection

Drawer rack type to facilitate  
the manipulation of the material



Insulated hinged lid with  
removable magnetic rubber  
gasket to guarantee a  
hermetic closure

Control Panel, microprocessor  
operating for all functions,  
acoustic and visual LED display

CFC-free refrigeration system

Compact sizes for better use of  
the laboratory space

### WHAT TO ORDER

	Code	Capacity (litres)	Dimensions	Weight (kg)	N° of doors	N° of shelves
Pannel	LED					
	<b>89503</b>	84	500x610x850 h mm	35	1	3 drawers
	<b>89505</b>	175	600x650x1400 h mm	50	1	5 drawers
	<b>89498*</b>	250	600x650x1850 h mm	90	1	8 drawers