

# Daily-Head

DISPOSABLE ASPIRATING HEAD FOR MICROBIOLOGICAL AIR MONITORING



- Disposable, triple packed, aspirating head
- Made in antistatic resin plastic
- Irradiated triple-wrapped bag for clean room
- Certification of sterilisation and conformity in five languages (English, Italian, French, German, Spanish)
- Suitable for "SAS Super", "DUO-SAS-360", "SAS-Isolator" air samplers.
- Suitable for any type of diameter 55 mm Contact Plates (RODAC).
- Ideal in combination with "Atrium" plate chamber for use in Clean Room and Isolator.
- Five years shelf-life.

86996



SAS Super



DUO SAS 360



SAS Isolator

**WHAT TO ORDER**

<b>86996</b>	"Daily-Head" disposable aspirating head - for 55 mm Contact Plates	50/box
<b>91214</b>	"Daily-Head" disposable aspirating head - for 50 mm Petri Dishes	50/box