

SAS DAILY-HEAD

DISPOSABLE ASPIRATING HEAD FOR SAS AIR MONITORING

- All analytical material that are used for microbiological testing should be sterile and, according to the cGLP, a document of certification should prove that it has been correctly sterilized.
- This rule is also valid for the aspiration head that is used with microbiological air samplers. Each cycle of air sampling should use a documented sterile aspirating head.
- The "Daily-Head", is a disposable aspirating head for microbiological air monitoring, which has been developed and produced by international pbi to simplify the workload of the microbiologist.
- The "Daily-Head", made in antistatic resin plastic, is double packed, irradiated and comes with certification of sterilisation, conformity and endotoxin level (LAL test)
- It has a shelf-life of five years and it is suitable for any type of SAS Super air sampler that uses diameter 55 mm Contact Plates (RODAC).



DAILY-HEAD - 86996

WHAT TO ORDER

CODE	DESCRIPTION
86996	"Daily-Head" for 55 mm Contact Plates Certified sterile - 50 pieces per box
91214	"Daily-Head" for 90 mm Petri Dishes - 50 pieces per box