

## Press Release

# STERILE CERTIFIED HEAD OF AIR SAMPLING

*Ref. 11/09 En*

**Milan, 25/03/2009**

The "Good Laboratory Practice" requests that all the used instruments for sampling and/or analytical tests are certified to guarantee a final correct analytical result. Each re-usable stainless steel or aluminium [aspirating head](#) of the "[SAS](#)" microbiological air samplers from "International PBI" is individually identified by an engraved serial number and complete of certificate of performance.



The "[Dispo-head](#)" plastic disposable aspirating head of the "SAS" is single packed, sterilised by gamma irradiation, and complete of lot number and certificate of irradiation. The certification should be part of the documentation to be presented to the FDA inspectors or local authorities during their visits.