

Press Release

MICROBIOLOGICAL AIR SAMPLER "SAS MULTI ISOLATOR"

Ref. 35/09 En

Milan, 23/10/2009

"[SAS Multi Isolator](#)" microbiological air sampler produced by International PBI in Milan (Italy) was developed specifically for Clean Room and Isolators. It consists of a programmable unit which remains outside the controlled area and different sampling head which uses 55 mm contact plates (RODAC) or Standard 90 mm Petri dishes. Ideal for critical air monitoring in isolator and barrier environments. The performances are the same as the [SAS Super 100](#) or [SAS Super 180](#). No disturbance of the unidirectional air flow and fast sampling (up to 1000 litres (= 1 cubic metre of air in six minutes). It is fully programmable from outside and is compatible with the most commonly used sanitizing and sterilizing agents.

